## Lipase Color

**CATALOGUE NUMBER:** 905-B **SIZE:** Reagent 5 x 30 mL  
Solvent 1 x 200 mL  
Activator 1 x 60 mL  
Calibrator 2 x 3 mL  
905-E Reagent 33 x 30 mL  
Solvent 1 x 1 L  
Activator 1 x 1 L  
Calibrator 6 x 3 mL

### INTENDED USE

For the quantitative measurement of pancreatic lipase activity in serum or plasma.

### TEST SUMMARY

Pancreatic lipase in serum or plasma is closely associated with pancreatic diseases. The activity of this enzyme is measured as an important marker for diagnosing pancreatic diseases and the associated monitoring of therapeutic effects. Pancreatic lipase measurement has been reported using titrimetric, turbidimetric, fluorometric, and colorimetric methodologies.<sup>1-6</sup>

The Pancreatic Lipase Color test is a colorimetric, kinetic assay that uses a clear substrate solution of 1,2-diglyceride that is a “natural” substrate. The assay is highly sensitive and specific for pancreatic lipase, using colipase and deoxycholate as activators.<sup>5</sup>

The assay shows excellent reproducibility and stability. Furthermore, the simplicity of this procedure makes it readily adaptable for use on automated analyzers.

### TEST PRINCIPLE

Serum lipase acts on a natural substrate, 1,2-diglyceride, to liberate 2-monoglyceride. This is hydrolyzed by monoglyceride lipase (a highly specific enzyme for monoglyceride) into glycerol and free fatty acid. Glycerol kinase acts on glycerol to form glycerol-3-phosphate, which is in turn acted on by glycerol-3-phosphate oxidase to generate hydrogen peroxide. Peroxidase converts the hydrogen peroxide, 4-Aminoantipyrine and TOOS (N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine) into a quinine dye. The rate of formation of the dye, measured as an increase in absorbance at 550 nm, is proportional to the lipase concentration in the sample.<sup>5</sup>

### REAGENTS

#### Composition

Component	Ingredients	Concentration
Lipase Color Reagent	1,2-Diglyceride (Egg)	1.1 mM
	Monoglyceride Lipase (Bacillus sp.)	0.88 U/mL
	Glycerol Kinase (S. canus)	<1.34 U/mL
	Glycerol-3-Phosphate Oxidase (Streptococcus sp.)	<40 U/mL
	TOOS	0.07%
	ATP (Bacterial)	0.66 mM
	Peroxidase (Horseradish)	<1.34 U/mL
	Colipase (Porcine)	<40 U/mL
	Buffer	pH 6.8
	Human Serum Albumin (Human)	0.27%
	Ascorbate Oxidase (Cucumber, Zucchini)	<2.66 U/L
	Stabilizers	

#### Composition (continued)

Component	Ingredients	Concentration
Lipase Color Solvent	Cholic Acid (Ox or Sheep)	5.3 mM
	Buffer	pH 6.8
	Sodium Azide	0.05%
Lipase Color Activator	Deoxycholate (Ox or Sheep)	36 mM
	4-Aminoantipyrine	0.12%
	Buffer	pH 8.7
	Sodium Azide	0.05%

The Lipase Color Calibrator contains human pancreatic lipase, serum albumin (bovine), and preservative.

### WARNINGS AND PRECAUTIONS FOR USE

Contains: Oxidase, glycerol phosphate

#### Danger

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

#### Prevention

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P284 Wear respiratory protection.

#### Response

P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P342 + P311 If experiencing respiratory symptoms: Call a POISON CENTRE/doctor.

#### Disposal

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- For *In Vitro* Diagnostic Use.
- Do not pipette by mouth.
- Do not use the calibrator after the expiration date printed on the label.
- Warning:** Human source material. Treat as potentially infectious. Each donor unit used in the preparation of this product has been tested by an FDA approved method and found non-reactive for HBsAg, HCV, HIV 1/2 and HIV-1 antigen. Because no known test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good laboratory practices using appropriate precautions.<sup>7</sup>
- The enzymes in triglyceride and cholesterol reagents may contaminate the Lipase Color Assay. To avoid contamination, ensure probes, cuvettes or tubes of automated analyzers are thoroughly washed between triglyceride or cholesterol assays and use of the Lipase Color Assay.
- Caution:** Contains sodium azide, which may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush drain with a large volume of water to prevent buildup. Dispose of in accordance with local, state, and federal regulations.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

- Lipase Color Reagent: Reconstitute the Reagent vial with the 30 mL Solvent, using a Class A volumetric pipette. Allow to stand a minimum of ten minutes at room temperature. Mix gently by inversion before use.
- Lipase Color Solvent: Supplied ready to use.
- Lipase Color Activator: Supplied ready to use.

- Lipase Calibrator: Reconstitute the calibrator vial with 3.0 mL of distilled, deionized, Type II water or equivalent using a Class A volumetric pipette<sup>8</sup>. Replace the cap/stopper and mix gently by inversion. Allow to stand for a minimum of 10 minutes at room temperature. Mix gently by inversion before use.

**DO NOT SHAKE**

Unopened Lipase Color Reagent, Solvent, Activator and Calibrator are stable until the expiration date shown on the label when stored at 2-8°C.

**DO NOT FREEZE**

Once reconstituted, the Lipase Color Reagent is stable for 28 days at 2-8°C.

**DO NOT FREEZE**

Once opened the Lipase Color Solvent and Activator are stable until the expiration date on the label when stored capped at 2-8°C.

**DO NOT FREEZE**

Once reconstituted, the Lipase Calibrator is stable up to 28 days at 2-8°C. Reconstituted stability of the calibrator may be extended by freezing the reconstituted calibrator at -20°C for up to 4 months. Freezing in the original vial is recommended.

**FREEZE ONLY ONCE**

**REAGENT DETERIORATION**

The following would indicate deterioration:

Presence of extreme turbidity.

**DISPOSAL**

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

**SPECIMEN**

Serum, EDTA-treated plasma or lithium and sodium heparinized plasma are the recommended sample types. Specimens should be collected as per the National Committee for Clinical Laboratory Standards Guideline H4-A3.<sup>9</sup>

**Serum:** Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection.<sup>9</sup> (Within 3 hours)

**Plasma:** Specimens may be collected in EDTA or lithium and sodium heparin. Centrifuge and remove the plasma as soon as possible after collection.<sup>9</sup> (Within 3 hours)

If the assay is not performed immediately, the serum and plasma must be refrigerated or frozen until use. If not analyzed promptly, specimens may be stored at 2-8°C for up to 3 weeks. If specimens need to be stored for more than 3 weeks, they may be preserved at -20°C or below for up to 3 months. Samples may be frozen once. Refer to NCCLS Document H18-A for further instructions on specimen collection, handling and storage.

**ANALYTICAL SPECIFICITY**

No significant interference was detected in the Lipase Color assay up to and including the concentrations stated below:

<b><u>Substance Tested</u></b>	<b><u>Concentration with no significant (±10%) interference</u></b>
Unconj. Bilirubin	20 mg/dL
Conj. Bilirubin	25 mg/dL
Hemoglobin	2000 mg/dL

Triglyceride	1000 mg/dL
Liposyn	1%
Glycerol	250 mg/dL
Ascorbic Acid	50 mg/dL

Samples containing interfering substances greater than the levels listed exhibited bias >10%. Dilutions are not recommended to reduce interferences.

- Refer to the work of Young et al<sup>10</sup> for a review of the effects of drugs on clinical laboratory tests.
- The enzymes in triglyceride and cholesterol reagents may contaminate the Lipase Color Assay. To avoid contamination, ensure probes, cuvettes or tubes of automated analyzers are thoroughly washed between triglyceride or cholesterol assays and use of the Lipase Color Assay.

**ANALYTICAL PROCEDURE**

**MATERIALS PROVIDED**

The Lipase Color Reagent, Solvent, Activator and Calibrator are required for the measurement of Lipase. The Lipase Color reagents are packaged and sold together. A configuration of each of the following items will be included in the package you receive.

<b><u>Description</u></b>	<b><u>Configuration</u></b>	<b><u>Catalog Number</u></b>
Lipase Color	Reagent 5 x 30 mL Solvent 1 x 200 mL Activator 1 x 60 mL Calibrator 2 x 3 mL	905-B
Lipase Color	Reagent 33 x 30 mL Solvent 1 x 1 L Activator 1 x 1 L Calibrator 6 x 3 mL	905-E

**MATERIALS REQUIRED (BUT NOT PROVIDED)**

- Spectrophotometer or other instrument capable of reading at 550 nm.
- Class A volumetric pipettes.
- Distilled or deionized water meeting specifications equivalent to USP purified water.
- Lipase control sera or quality control material (See "Quality Control").

**TEST CONDITION**

All analyzer applications should be validated in accordance with CLIA recommendations. For assistance with applications on automated analyzers, please contact Sekisui Diagnostics Technical Services at 800-565-0265.

**CALIBRATION**

The Lipase Color Calibrator is required for calibration. Refer to the instrument manual for analyzer specific calibrator procedures and for guidance in determining calibration frequency. The calibrator value can be found on the calibrator vial label.

**QUALITY CONTROL**

Reliability of test results should be routinely monitored with quality control materials or serum pools that reasonably represent performance on patient specimens. Controls or serum pools should be run with each assay to ensure that the reagents are functioning properly and that correct procedures have been followed. Quality control materials are intended for use only as monitors of accuracy and precision. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range confirm procedures were performed correctly and follow normal troubleshooting measures. If the problem persists call Sekisui Diagnostics Technical Services at 800-565-0265.

Quality control requirements should be established in accordance with local, state and/or Federal regulations or accreditation requirements.

## REFERENCE INTERVALS

A normal range study was performed using the Lipase Color assay. A serum range of 21-67 U/L (138 healthy donors) was obtained on the Roche/Hitachi 911 system. Ranges were calculated as recommended by NCCLS guideline C28-A.<sup>11</sup> These results were obtained using a specific lot of reagent. It is recommended that each laboratory establish the normal range for its patient population.

## PERFORMANCE CHARACTERISTICS

### RESULTS

1. Unit Definition: One unit (U) is defined as the amount of enzyme activity which liberates 1  $\mu$ mole of 2-monoacylglyceride from 1,2-diglyceride per minute at 37°C. Lipase activity is expressed in U/L.
2. Patient results may be reported in U/L. To convert from conventional units to S.I. units, multiply the conventional units by  $1.65 \times 10^{-8}$  Katal/U.<sup>5</sup> U/L  $\times (1.6 \times 10^{-8}$  Katal/U) = Katal/L Lipase.

### REPORTABLE RANGE

#### Limit of Detection

Limit of detection of the Lipase Color assay, quantified as 2 S.D. plus the mean of twenty replicate measurements of saline, is 2 U/L on a Roche/Hitachi 911 analyzer.

#### Linearity

The Lipase Color method is linear to 750 U/L lipase on the Roche/Hitachi 911 analyzer. If the result is greater than 750 U/L, dilute with saline, multiply the result by the dilution factor to obtain the lipase activity of the sample.

### PRECISION STUDIES

Within-run precision of the Lipase Color assay was determined using three frozen spiked human lipase pools. Each run consisted of 20 replicate samples. Within-run precision studies produced the following results on the Roche/Hitachi 911 Analyzer.

Serum Pool	Low	Mid	High
n	20	20	20
Mean (U/L)	33	118	269
S.D. (U/L)	0.8	1.5	2.1
CV (%)	2.4	1.2	0.8

Between run precision of the Lipase Color assay was determined using three frozen spiked human lipase pools tested in duplicate, once or twice per day, for 15 days.

Serum Pool	Low	Mid	High
n	20	20	20
Mean (U/L)	34	120	275
S.D. (U/L)	1.5	2.7	6.3
CV (%)	4.4	2.3	2.3

## ACCURACY

Accuracy of the Lipase Color method, performed using the Roche/Hitachi 911 Analyzer was verified by comparison to the Dade Lipase method, performed on a Dade Dimension analyzer, producing the following results:

Method	Sekisui Diagnostics	Dade Lipase Assay
Lipase Assay	n	40
Mean (U/L)	54	260
Regression Analysis	$y=0.44(x)-62.07$	
Correlation Coefficient	$r=0.97$	
Reference Range	21-67 U/L	114-286 U/L

All trademarks, brands, product names and trade names are the property of their respective companies.

The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

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**IVD**

*In vitro* diagnostic medical device  
Appareil médical de diagnostic *in vitro*.

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Caution, consult accompanying document  
Prudence, consulter les documents joints

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1. Imamura S., Misaki H., "A sensitive method for assay of lipase activity by coupling with  $\beta$ -oxidation enzymes of fatty acid." Selected Topics in Clinical Enzymology; 2:73 (1984).
2. Imamura S., et al, "A sensitive method for assay of Lipase activity using 1,2-diglyceride as substrate and coupling with  $\beta$ -oxidation enzymes of fatty acid as an indicator reaction." Collection of summaries of lectures in the 126th general meeting of Kinki branch, analytical section, Japan Society of Clinical Chemistry; p.11-31 (1986).
3. Hayashi C., et al, "Assay methods for human lipases." Clinical Examination, Instrument and Reagent, 2:225 (1986).
4. Kitaura S., et al, "Properties of Monoglyceride Lipase Produced by Thermophillic Bacteria." Seikagaku 1988 Aug; 60 (8): 848.
5. Imamura S., et al, "An Enzymatic Method Using 1,2-Diglyceride for Pancreatic Lipase Test in Serum." Clin. Chem. 1989; 35 (6): 1126.
6. Tietz NW, "Clinical Guide To Laboratory Tests", 3rd ed., Philadelphia, PA; WB Saunders Co.; 364 (1995).

7. Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 1999.
8. National Committee for Clinical Laboratory Standards, Preparation and Testing of Reagent Water in the Clinical Laboratory-Third Edition: approved Guideline NCCLS Document C3-A3, 1997.
9. National Committee for Clinical Laboratory Standards, "Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture", Approved Standard, Third Edition, NCCLS publication H4-A3, Villanova, PA (1991).
10. Young DS, Effects of Drugs on Clinical Laboratory Tests, 4th Edition, AACC Press, Washington, DC; 3-398 to 3-400 (1995).
11. National Committee for Clinical Laboratory Standards, "How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline", NCCLS Document C28-A. Vol. 15, No.4, June 1995.

©2011 Sekisui Diagnostics, LLC - All rights reserved.

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©2011 Sekisui Diagnostics, LLC - Tous droits réservés.