

MATERIAL SAFETY DATA SHEET

I PRODUCT IDENTIFICATION AND USE

SUPPLIERS NAME: Diagnostic Chemicals Limited
16 McCarville St
West Royalty Industrial Park
Charlottetown, P.E.I.
Canada C1E-2A6
902-566-1396

CHEMICAL/PRODUCT NAME: **Direct HDL Calibrator**

CATALOGUE NUMBER: **SE-248**

NOTE: Physical and health hazard information on reagent mixtures has not been determined. Any physical and health information noted is based 1) on evaluation of data for the pure ingredients, and 2) concentration of ingredients as packaged.

PRODUCT DESCRIPTION: To calibrate in-vitro diagnostic test for HDL cholesterol

II HAZARDOUS INGREDIENTS

HDL Cholesterol Calibrator: n.ap. ACGIH TLV (NE = Not Established) CAS

Human sourced material. Although the product was found negative to hepatitis B surface antigen and antibodies to HIV 1&2 and HCV, it should be handled at Biosafety Level 2 as recommended for any potentially infectious serum and specimen.

III TOXICOLOGICAL PROPERTIES

ROUTE OF ENTRY/EXPOSURE

SKIN CONTACT: [x] EYE CONTACT: [x] INHALATION: [x] INGESTION: [x]
[x] SKIN ABSORPTION [x]

EFFECTS OF ACUTE EXPOSURE:

SKIN CONTACT: May cause irritation.

EYE CONTACT: May cause irritation.

INGESTION: May be harmful if ingested.

INHALATION: May cause irritation to mucous membranes and upper respiratory tract.

EFFECTS OF CHRONIC EXPOSURE: n.ap.

EVIDENCE OF CARCINOGENICITY, TERATOGENICITY, REPRODUCTIVE TOXICITY, MUTAGENICITY: n.av.

SENSITIZATION TO PRODUCT: n.av.

SYNERGISTIC PRODUCTS: n.av.

LD50: n.av.

IV EMERGENCY AND FIRST AID MEASURES

SKIN: In case of contact, immediately flush area with plenty of water.

EYES: Immediately flush eyes with plenty of water for at least 15 minutes and get medical attention.

ORAL: If ingested drink plenty of fluids and contact a physician.

INHALATION: If inhaled remove victim to fresh air. If not breathing give artificial respiration. If breathing is difficult give oxygen, and call a physician.

V PHYSICAL DATA

PHYSICAL STATE: Solid

APPEARANCE AND ODOUR: Pale yellow powder

ODOUR THRESHOLD: n.av.

VAPOUR PRESSURE: n.av.

VAPOUR DENSITY: n.av.

SPECIFIC GRAVITY: n.av.

EVAPORATION RATE: n.av.

MELTING POINT: n.av.

BOILING POINT: n.av.

SOLUBILITY IN WATER: Soluble

pH: n.av.

COEFFICIENT OF OIL/WATER DISTRIBUTION: n.av.

VI FIRE AND EXPLOSION HAZARD DATA

FLAMMABLE: No

EXTINGUISHING MEDIA: As appropriate for fire in surrounding materials

FLASH POINT: n.ap.

UPPER FLAMMABLE LIMIT: n.ap.

LOWER FLAMMABLE LIMIT: n.ap.

AUTOIGNITION TEMPERATURE: n.ap.

HAZARDOUS COMBUSTION PRODUCTS: n.av.

EXPLOSION DATA:

SENSITIVITY TO MECHANICAL IMPACT: n.ap.

SENSITIVITY TO STATIC DISCHARGE: n.ap.

UNUSUAL FIRE AND EXPLOSION HAZARDS: Unknown

VII REACTIVITY DATA

STABLE: Yes

HAZARDOUS POLYMERIZATION: Will not occur

INCOMPATIBILITY WITH OTHER SUBSTANCES: n.av.

HAZARDOUS DECOMPOSITION PRODUCTS: n.av.

VIII PREVENTIVE MEASURES

LOCAL EXHAUST: Not required
PROTECTIVE CLOTHING: Lab coat
PROTECTIVE GLOVES: Rubber
EYE PROTECTION: Safety glasses

OTHER PRECAUTIONS: Avoid contact and inhalation, do not get in eyes, on skin or clothing. Do not pipet by mouth. Wash contaminated clothing before reuse.

RESPIRATORY PROTECTION: None normally required.

LEAK AND SPILL PROCEDURE: Absorb small leaks or spills with sponge, mop up large spills with plenty of soap and water.

WASTE DISPOSAL: Dilute with a large volume of water and dispose of into sewer system, in accordance with local regulations.

STORAGE REQUIREMENTS: 2 -8°C

SHIPPING INFORMATION: Not regulated

IX PREPARATION OF THE M.S.D.S.

PREPARED BY: Denise Bustard
PHONE NO: 902-566-1396
DATE: February 16, 2001

** n.av. not available n.ap. not applicable*

The information provided on this Material Safety Data Sheet is provided in the interest of promoting safe handling of the material. While this information is believed to be correct DCL makes no warranty in respect to any of the information disclosed. Observe all federal, provincial, state and local laws concerning health and pollution.

MSDSSE248-1